

# Recording, analysis and visualization of data from Advacam miniPIX (EDU) silicon pixel detector

Günter Quast, February 2026

## mPIXdaq Data acquisition and analysis for *miniPIX (EDU)* pixel detector

Vers. 1.0.1, February 2026

The miniPIX EDU is a camera for radiation based on the Timepix pixel read-out chip with 256x256 radiation-sensitive pixels of 55x55 $\mu\text{m}^2$  area and 300 $\mu\text{m}$  depth each. The chip is covered by a very thin foil to permit  $\alpha$  and  $\beta$  radiation to reach the pixels. The device is enclosed in an aluminum housing with a USB 2.0 interface. The sensor chip is covered with a thin foil and is very fragile; this area should be protected with a cover if not measuring  $\alpha$  radiation.

The device provides two-dimensional images of particle traces in the sensitive detector material. The high spatial resolution compared to the typical range of particles in silicon is useful to distinguish the different types of radiation and to measure their deposited energies.  $\alpha$ -particles are completely absorbed and deposit all of their energy in the sensitive area, allowing usage of the device as an energy spectrometer.

The vendor provides a ready-to-use program for different computer platforms as well as a software-development kit for own applications.

The code provided here is a minimalist example to read out single frames, i.e. a full set of 256x256 pixel energies accumulated over a given, fixed time interval. Each frame is displayed as an image with a logarithmic color scale representing the deposited energy in each pixel.

The analysis of the recorded signals, i.e. clustering of pixels, energy determination and visualization, is achieved with standard open-source tools for data analysis. It is therefore well-suited to provide high-school or university students with detailed insights into principles of the interaction of radiation with matter and to enable them to carry out their own studies.

## Getting ready for data taking

This code has been tested on *Ubuntu*, *openSuse*, *Fedora*, on Windows 64bit with *Python3.7.9* and on *Raspberry Pi* for the 32- and 64-bit versions of *OS12* and *OS13*. Other Linux distributions should not pose any unsurmountable problems.

On MS Windows, the libraries provided by the vendor only support *Python* vers. 3.7.9; such a rather historic version can be set up using e.g. the *miniconda* framework.

The code also supports devices other than the miniPIX EDU if the configuration files are available and copied to the *factory/* directory in the *pypixet Python* interface.

To get started, follow the steps below:

- Get the code from github  

```
git clone https://github.com/GuenterQuast/mPIXdaq.
```

This repository includes the *Python* code and a minimalistic set of libraries provided by Advacam.
- Next, `cd` to the `mPIXdaq` directory you just downloaded.
- Set up the USB interface of your computer to recognize the miniPIX detector:  

```
sudo install_driver_rules.sh
```

 (to be done only once), then connect the *miniPIX* to your computer.

The package may also be installed in your virtual python environment:

- ```
python -m pip install .
```

Now everything is set up to enjoy your *miniPIX*. Just run the *Python* program from any working directory by typing

```
run_mPIXdaq.py.
```

If you plan to record data, note that the path to the output file is relative to the current working directory.

*Note* also that the *pypixet* initialization is set up to write log-files and configuration data to the directory */tmp/mPIX/*.

It is also worth mentioning that on some systems the current directory, ".", needs to be contained in the LD\_LIBRARY\_PATH so that the Advacam *Python* interface *pypixet* finds all its C libraries and configuration files. This is also done in the *Python* script *run\_mPIXdaq.py* by temporarily modifying the environment variable LD\_LIBRARY\_PATH if necessary and then restarting to execute the *Python* code in the new environment.

## Running the example script

Available options of the *Python* example to steer data taking and data archival to disk are shown by typing

*run\_mPIXdaq.py --help*, resulting in the following output:

```
usage: run_mPIXdaq.py [-h] [-v VERBOSITY] [-o OVERLAY] [-a ACQ_TIME]
                    [-c ACQ_COUNT] [-f FILE] [-w WRITEFILE] [-t TIME]
                    [--circularity_cut CIRCULARITY_CUT]
                    [--flatness_cut FLATNESS_CUT] [-r READFILE]
```

read, analyze and display data from miniPIX device

options:

```
-h, --help                show this help message and exit
-v VERBOSITY, --verbosity VERBOSITY
                           verbosity level (1)
-o OVERLAY, --overlay OVERLAY
                           number of frames to overlay in graph (10)
-a ACQ_TIME, --acq_time ACQ_TIME
                           acquisition time/frame in seconds (0.5)
-c ACQ_COUNT, --acq_count ACQ_COUNT
                           frame count for readout via callback (20)
-f FILE, --file FILE      file to store frame data
-w WRITEFILE, --writefile WRITEFILE
                           file to write cluster data
-t TIME, --time TIME      run time in seconds (36000)
--circularity_cut CIRCULARITY_CUT
                           cut on circularity for alpha detection (0.5)
--flatness_cut FLATNESS_CUT
                           cut on flatness for alpha detection (0.6)
-p PRESCALE, --prescale PRESCALE
                           prescaling factor for frame analysis
-r READFILE, --readfile READFILE
                           file to read frame data
-b BADPIXELS, --badpixels BADPIXELS
                           file with bad pixels to mask
```

The default values are adjusted to situations with low rates, where frames from the *miniPIX* with an exposure time of *acq\_time* = 0.5 s are read. For the graphics display, *overlay* = 10 recent frames are overlaid, leading to a total integration time of 5 s. These images represent a two-dimensional pixel map with a color code indicating the energy measured in each pixel.

The miniPIX EDU version, in particular, may suffer from a large number of dead or noisy pixels, and therefore they should be masked by providing a file with the pixel indices to be ignored. The default file name is *snxxx\_badpixels.txt* in the working directory, where *snxxx* is the serial number of the miniPIX device. Alternatively a file name may be specified using the *-b* or *--badpixels* option.

Collected frame data may be directly written to disk, if a filename is given using the `-f` or `--file` option. Two formats are foreseen at present, storage of the two-dimensional frames as *numpy*-arrays (file extension `.npy`) or as lists of pixel indices and energy values in *yaml*-format (file extension `.yaml`). If no suffix for the filename is given, the default behavior is writing a `.yaml` file. The `.yaml` format also permits storing meta data like parameters of the data acquisition, the properties of the sensor and the list of bad pixels. To save space, the output files may be compressed with *zip* or *gzip*.

The same formats are recognized when reading back files using the `-r` resp. `--readfile` options.

In addition, files written with the *Pixet* program of Advacam in `.clog` format can be used as an input to this software package.

Data analysis consists of clustering of pixels in each frame and determination of cluster parameters, like the number of pixels, energy of clusters, and the shapes of the cluster areas and of the energy distribution over the pixels in the clusters.

The shape of the cluster area is encoded in a quantity called circularity. For discrimination between linear and round clusters, a cut is controlled by the parameter `circularity_cut` ranging from 0. for perfectly linear to 1. for perfectly circular clusters. Technically, the covariance matrix of the cluster area is calculated, and the circularity is defined as the square root of the ratio of the smaller and the larger of the two eigenvalues of the covariance matrix. This simple procedure already provides a good separation of  $\alpha$  and  $\beta$  particles and of isolated pixels not assigned to clusters. The latter ones have a high probability of being produced in interactions of  $\gamma$  radiation (keV-photons), while electrons from  $\beta$  radiation or from  $\gamma$  interactions produce long traces.  $\alpha$  particles produce large, circular clusters due to their very high ionization loss in the detector material.

A further, very sensitive variable is the variance of the energy distribution in the clusters. For  $\alpha$  particles, this distribution peaks at the centre and steeply falls off towards the boundary, leading to a very small variance. A small ratio of the variances of the energy distribution and of the area covered by pixels is therefore a very prominent signature of  $\alpha$  particles. The cut separating peaking and flat signatures is controlled by the parameter `flatness` with values between 0 and 1.

Properties of clusters, including a list of contributing pixels and their energy values, are optionally written to a file in *yaml* format (file extension `.yaml`) for later off-line analysis. A more compact version in `.csv` format, containing just the cluster properties, is also available.

A *Jupyter* notebook, *analyze\_mPIXclusters.ipynb*, illustrates an example analysis based on such file formats.

To test the software without access to a *miniPIX* device or without a radioactive source, a file with recorded data is provided. Use the option `--readfile data/BlackForestStone.yaml.gz` to start a demonstration. Note that the analysis of the recorded pixel frames is done in real time and may take some time on slow computers.

## Parameter settings for data acquisition

The optimal choice of parameters, in particular the values of exposure time, overlay of frames for the graphical display depend very much on the use case.

In scenarios with low rates aiming at a demonstration of the capabilities of a modern particle detector like the *miniPIX* for outreach purposes, it is most useful to mimic the behavior of a cloud chamber, i.e. particle traces appear, remain visible for some time and then disappear again. This can be achieved by setting an acquisition time of 0.2 s and an overlay of 10 frames. Particles tracks then remain visible on screen for 2 s. The command to run in this mode is:

```
run_mPIXdaq -a 0.2 -o10
```

If the goal is to efficiently record particle tracks with low read-out dead-time, it is not necessary to optimize the visual impression, because the graphical display is only used to ensure the quality of the recorded data. To reduce processing overheads, the fraction of events being analyzed and displayed may be pre-scaled. In the example below, only every fourth frame is analyzed and displayed:

```
run_mPIXdaq -a 0.2 -p4 -o1
```

In high-rate scenarios exceeding 100 particles/s the read-out efficiency for the *miniPIX* becomes a concern. The USB 2 transfer and initialization overheads take about 25 ms per frame. In practice, this means that about 20 frames/s of 25 ms exposure time each can be transferred at a read-out dead-time of 50%. The number of objects per frame should not exceed about 100 in order to avoid overlaps between clusters. So, in practice, signatures of 2000 particles/s can be handled, which is clearly sufficient for most laboratory experiments with rather weak radioactive sources that

comply with radiation protection regulations. Read-out of the miniPIX is fastest in callback mode, when the driver is initialized to call a function for data retrieval whenever a new frame is ready to be transferred. To receive *acq\_count* frames, only one initialization step is necessary. The exposure time of each frame is given by the value of *acq\_time*. To achieve maximum read-out speed for such a high-rate scenarios, start data-acquisition with the command:

```
run_mPIXdaq -a 0.025 -c50 -p10 -o1
```

With these settings, only one tenth of the frames is analyzed and displayed; the recorded frame rate is 20 Hz, while the read-out dead-time indeed turns out to be 50% (measured on a Raspberry Pi 5).

If higher higher read-out rates up to the nominal 40 frames/sec are needed, data can be recorded using the *Pixet* (basic) program delivered by Advacam together with the device. Select tracking mode, the required exposure-time per frame and the number of frames from the graphical interface, then press the record button, and after completion save the acquired frames using the *save* button selecting *.clog* as the file format. This format is understood by *mPIXdaq* and can be read in via the *-readfile* option,

thus allowing you to perform cluster analysis and writing files in the same formats as implemented in *mPIXdaq*. In case your miniPIX device suffers from noisy pixels, a bad-pixel map can be specified using the *-badpixel* option.

## Implementation Details

The default data acquisition is based on the function *doSimpleAcquisition()* from the *Advacam Python* API in callback mode, i.e. *acq\_counts* frames with an adjustable accumulation time *acq\_time* are read from the miniPIX device successively.

The chosen readout mode is “ToT” (“time over threshold”, *PX\_TPXMODE\_TOT*). This quantity shows good proportionality to the deposited energy at high signal values, but exhibits a non-linear behavior for small signals near the detection threshold of the *miniPIX*. Calibration constants are stored on the miniPIX device for each pixel, which are used to provide deposited energies per pixel in units of keV.

The relevant libraries for device control are provided in directories *advacam\_<arch>* for x86\_64 Linux, arm32 and arm64 and for Macintosh arm64 and MS Windows architectures. The contents of a typical directory is:

```
__init__.py    # package initialization
pypixet.so     # the Pixet Python interface
minipix.so     # C library for pypixet
pxcore.so      # C library for pypixet
pixet.ini      # initialization file, in same directory as pypixet
factory/       # initialization constants
```

Note that the copyright of these libraries is held by Advacam. The libraries may be downloaded from their web page, ADVACAM DWONLOADS. They are provided here as *Python* packages for some platforms for convenience.

## Data Analysis

The analysis shown in this example is intentionally very simple and based on standard libraries and functions. Clustering of pixels is performed by finding connected regions in the pixel image with *scipy.ndimage.label()*. The shape of the clusters is determined from the ratio of the smaller and the larger one of the two eigenvalues of the covariance matrix calculated from the *x* and *y* coordinates of the pixels in a cluster using *numpy.cov()*. For circular clusters, as typically produced by  $\alpha$  radiation, this ratio is close to one, while it is almost zero for the longer traces from  $\beta$  radiation. In addition, the shape of the energy distribution is considered, which shows a sharp maximum at the center for  $\alpha$  particles but is rather flat otherwise.

The figure below shows the graphical display with a pixel image and the typical distributions of the pixel and cluster energies and the number of pixels per cluster. The source used was a weakly radioactive stone from the Black Forest containing a small amount of Uranium and its decay products. The pixel map shown in the figure was sampled over a time of five seconds. The histogram in the lower-right corner demonstrates that the cluster types of different types of radiation are well separated:  $\alpha$  rays in the green band with relatively low numbers of pixels per cluster, electrons ( $\beta$ ) as long tracks with large numbers of pixels per cluster and rather low energies. Single pixels not associated to clusters mostly originate from  $\gamma$  rays. Some of the electron tracks with typically low energies also stem from photon interactions in the detector material (via the Compton process).




Abbildung 1: The graphical display of miniPIXdaq

The frame collection time is chosen to be on the order of seconds, so that analysis results can be displayed in real-time on a sufficiently fast computer including the Raspberry Pi 5. This is suitable for investigations of natural radiation as emitted by minerals like Pitchblend (=Uraninit), Columbit, Thorianit and others. Radon accumulated from the air in basement rooms on the surface of an electrostatically charged ballon also work fine.

For applications at higher rates, the analysis may have to be done off-line by reading data from recorded files. The option `-prescale` can be used to limit frame and cluster analysis to a subset of the recorded frames while still allowing to record all data to file with sufficiently low dead time.

In a future version of this program an option to use multiple cores for the analysis task may be provided.

## Output files and formats

*miniPIXdaq* offers the possibility to write raw or clustered data for further off-line analysis, which is a valuable asset for use of the software in physics laboratory courses. This feature enables students to develop their own strategies for the analysis of data recorded in the student lab.

The default output-format is *yaml*, offering human readability, clear structure and modularity of data blocks. Output files are written sequentially as text files while data acquisition is progressing. Files may be compressed using *gzip* or *zip* to obtain more compact representations of the recorded data.

**Frame data** and **cluster data** including some meta data are stored in *yaml* structures with the keys *meta\_data:*, *deviceInfo:*, *bad\_pixels:* and *eor\_summary:*. Frame or cluster data are stored as lists of pairs of pixel indices and pixel energies for all pixels with non-zero energy under the keys *frame\_data:* or *cluster\_data:*, respectively. These files can be loaded into a *Python* dictionary using the python code

```
> yaml_dict = yaml.load(open('<filename>', 'r'), Loader=yaml.CLoader)
```

and frame or cluster data unpacked into a *Python* list via

```
> list_of_framedata = yaml_dict["frame_data"] or  
> list_of_clusterdata = yaml_dict["cluster_data"], respectively.
```

Frame data in text form or as *zipped* or *gzipped* files can be used as input to *mPIXdaq* via the `-r` or `-readfile` options.

*list\_of\_clusterdata* is a list of two lists, the first one containing cluster properties and the second one indices and energies of contributing pixels, i.e.

```
> list_of_clusterproperties[i] = yaml_dict["cluster_data"][i][0] and  
> list_of_clusterpixels[i] = yaml_dict["cluster_data"][i][1],
```

A *jupyter* notebook *analyze\_mPIXclusters.ipynb* is distributed as part of the package and illustrates how to read and interpret cluster data.

The keys of the variables in *list\_of\_clusterproperties* are

```
>['time', 'x_mean', 'y_mean', 'n_pix', 'energy', 'var_mx', 'var_mn', 'angle', 'xE_mean', 'yE_mean', 'varE_mx',  
'varE_mn']
```

It is also possible to store the cluster properties in simple *.csv* format by explicitly specifying the file extension:  
`run_mPIXdaq <options> -w <name>_clusters.csv.`

Histograms displayed in the on-line graphical display may be saved using the control buttons in the *matplotlib* window.

## Sensor Details

The miniPIX (EDU) is based on the *Timepix* hybrid silicon pixel device, consisting of a semiconductor detector chip segmented into 256 x 256 square pixels with a pitch of 55 mm that is bump-bonded to the readout chip. Each element of the pixel matrix is connected to its own preamplifier, discriminator and digital counter integrated on the readout chip.

The built-in *Medipix2* variant of the chip is operated in the so-called “frame mode”, i.e. all pixels are read out at the same time, providing one frame consisting of the deposited energies per pixel collected during the acquisition time. If operated in time-over-threshold (ToT) mode, returned pixel readings represent the time the signal is over a given threshold in counts of the chip clock (appr. 10 MHz). *ToT* is linearly related to the energy deposition for large deposits exceeding 50 keV. The functional dependence on the deposited energy  $E$ , including threshold effects, is approximated by the following function

$$ToT = a E + b - c/(E - t)$$

Approximate values of the calibrations constants are  $a = 1.6$ ,  $b=23$ ,  $c=23$  and  $t=4.3$ . Each pixel has its individual calibration stored on the chip, which is optionally applied to obtain pixel readings in units of keV. The calibration is reliable up to pixel energies of one MeV. Higher pixel energies may result when frames with short acquisition time are summed up. For details, see the article by J. Jakubek, *Precise energy calibration of pixel detector working in time-over-threshold mode*, NIM A 633 (2011), 5262-5265\*.

## Package Structure

This package consists of one *Python* file with several classes providing the base functionality. As mentioned above, it relies on Advacam libraries for setting-up and reading the sensor. Other dependencies are well-known libraries from the “Python” eco-system for data analysis:

- numpy
- matplotlib
- scipy.ndimage.label
- numpy.cov
- numpy.linalg.eig

The classes and scripts of the package are

- class miniPIXdaq
- class frameAnalyzer
- class miniPIXvis
- class runDAQ
- class bhist
- class scatterplot
- package script run\_mPIXdaq.py

Details on the interfaces are given below.

```
class miniPIXdaq:
    """Initialize and readout miniPIX device

    After initialization, the __call__() method of this class is executed
    in an infinite loop, storing data from the device in a ring buffer.
    The current buffer index is sent to the calling process via a Queue
    (dataQ, an instance of threading.Queue()). The loop ends when the flag
    endEvent (an instance of threading.Event() in class mpixControl) is set.

    Args:

        - ac_count: number of frames to read successively
        - ac_time: acquisition (=exposure) time per frame

    Queues for communication and synchronization

        - dataQ: Queue to transfer data
        - cmsQ: command Queue

    Data structure:

        - fBuffer: ring buffer with recent frame data

    """

class frameAnalyzer:
    """Analyze frame data and produce a list of cluster objects,

    Args:  2d frame data, as obtained from miniPIXdaq.__call__()
    """
```

Output:

```
pixel_clusters: a list of tuples of format
    (x, y), n_pix, energy, (var_mx, var_mn), angle, (xEm, yEm), (varE_mx, varE_mn)
with cluster properties
```

Helper functions to store analysis results are included as static methods

Another static method, `cluster_summary()` is particularly useful for on-line monitoring of incoming data and provides a summary of the properties of clusters in a pixel frame, returning

- `n_clusters`: number of multi-pixel clusters
- `n_cpixels`: number of pixels per cluster
- `circularity`: circularity per cluster (ranging from 0. for linear, 1. for circular)
- `flatness`: ratio of maximum variances of pixel and energy distributions in clusters
- `cluster_energies`: energy per cluster
- `single_energies`: energies in single pixels

```
class miniPIXvis:
```

```
    """Display of miniPIX frames and histograms for low-rate scenarios,
    where on-line analysis is possible and animated graphs are meaningful
```

```
    Animated graph of (overlaid) pixel images, number of clusters per frame
    and histograms of cluster properties
```

```
    Args:
```

- `npix`: number of pixels per axis (256)
- `nover`: number of frames to overlay
- `unit`: unit of energy measurement ("keV" or "µs ToT")
- `circ`: circularity of "round" clusters (0. - 1.)
- `flat`: flatness of energy distribution of pixels in clusters (0. - 1.)
- `acq_time`: accumulation time per read-out frame
- `prescale`: prescale factor for frame analysis

Objects of these classes are instantiated by the class `runDAQ`.

This class also accepts the command-line arguments to set various options, as already described above.

```
class runDAQ:
```

```
    """run miniPIX data acquisition, analysis and real-time graphics and data storage
```

```
class to handle:
```

- command-line arguments
  - initialization of miniPIX device of input file
  - real-time analysis of data frames
  - animated figures to show a live view of incoming data
  - event loop controlling data acquisition, data output to file and graphical display
- ```
"""
```

Two helper classes implement 1d and 2d histogram functionality for efficient and fast animation using methods from `matplotlib.pyplot`.

```
class bhist:
```

```
    """one-dimensional histogram for animation, based on bar graph
    supports multiple data classes as stacked histogram
```

```
    Args:
```

- \* `data`: tuple of arrays to be histogrammed



```

        * bindegges: array of bin edges
        * xlabel: label for x-axis
        * ylabel: label for y axis
        * yscale: "lin" or "log" scale
        * labels: labels for classes
        * colors: colors corresponding to labels
    """

class scatterplot:
    """two-dimensional scatter plot for animation, based on numpy.histogram2d.
    The code supports multiple classes of data and plots a '.' in the corresponding color in every non-zero

    Args:
        * data: tuple of pairs of coordinates (([x], [y]), ([], []), ...)
          per class to be shown
        * binedges: 2 arrays of bin edges ([bex], [bey])
        * xlabel: label for x-axis
        * ylabel: label for y axis
        * labels: labels for classes
        * colors: colors corresponding to labels
    """

```

A package script `run_mPIXdaq` is provided as an example to tie everything together into a running program. Because the ADVACAM *Python* interface (`pypixet.so`) expects C-libraries and configuration files in the very same directory as the Python interface `pypixet.so` itself, some tricky manipulation of the environment variable `LD_LIBRARY_PATH` is needed to ensure that all libraries are loaded and the *miniPIX* is correctly initialized.

```

#!/usr/bin/env python
#
# script run_mPIXdaq.py
# run mpixdaq example with data acquisition, on-line analysis and visualization
# of pixel frames and histogramming

import os, platform, sys

# on some Linux systems, pypixet requires '.' in LD_LIBRARY_PATH to find C-libraries
# - add current directory to LD-LIBRARY_PATH
# - and restart python script for changes to take effect

path_modified = False
if 'LD_LIBRARY_PATH' not in os.environ and platform.system() != 'Windows':
    os.environ['LD_LIBRARY_PATH'] = '.'
    path_modified = True
    print(" ! temporarily added '.' to LD_LIBRARY_PATH !")
    # restart script in modified environment
    try:
        os.execv(sys.argv[0], sys.argv)
    except Exception as e:
        sys.exit('!!! run_mPIXdaq: Failed to Execute under modified environment: ' + str(e))

# get current working directory (before importing minipix libraries)
wd = os.getcwd()

if os.name == 'nt':
    # special hack for windows python 3.7: load pypixet and DLLs
    import mpixdaq.advacam_win64.pypixet as pypixet

from mpixdaq import mpixdaq # this may change the working directory, depending on system

```

```
# start daq in working directory
rD = mpixdaq.runDAQ(wd)
rD()
```

It is also possible to start the script as a *Python* module:

```
python -m mpixdaq
```