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```
> run_mPIXdaq.py -h
usage: run_mPIXdaq.py [-h] [-v VERBOSITY] [-o OVERLAY] [-a ACQ_TIME] [-c ACQ_COUNT] [-f FILE] [-w WRITEFILE] [-t TIME]
                        [--circularity_cut CIRCULARITY_CUT] [--flatness_cut FLATNESS_CUT] [-p PRESCALE] [-r READFILE] [-b BADPIXELS]
                        [--no-kbdControl] [--no-guiControl]
```

Read, analyze & display miniPIX data

options:

```
-h, --help            show this help message and exit
-v VERBOSITY, --verbosity VERBOSITY
                        verbosity level (1)
-o OVERLAY, --overlay OVERLAY
                        number of frames to overlay in graph (4)
-a ACQ_TIME, --acq_time ACQ_TIME
                        acquisition time/frame in seconds (0.5)
-c ACQ_COUNT, --acq_count ACQ_COUNT
                        frame count for readout via callback (20)
-f FILE, --file FILE  file to store frame data
-w WRITEFILE, --writefile WRITEFILE
                        file to write cluster data
-t TIME, --time TIME  run time in seconds (36000)
--circularity_cut CIRCULARITY_CUT
                        cut on circularity for alpha detection (0.5)
--flatness_cut FLATNESS_CUT
                        cut on flatness for alpha detection (0.6)
-p PRESCALE, --prescale PRESCALE
                        prescaling factor for frame analysis
-r READFILE, --readfile READFILE
                        file to read frame data
-b BADPIXELS, --badpixels BADPIXELS
                        file with bad pixels to mask
--no-kbdControl        switch off keyboard control
--no-guiControl        switch off gui control
```

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```
*==* Contents of file ../data/Radon_clusters.yaml.gz  
Data written on Sun Apr 12 23:36:34 2026 with device EduMiniPIX TSWB sn: 2987  
16313 frames of 0.5s exposure time  
47227 clusters -> rate = 5.79 Hz
```

```
*==* cluster features:  
Index(['time', 'x_mean', 'y_mean', 'n_pix', 'energy', 'e_mx', 'x_mn', 'y_mn',  
      'w', 'h', 'var_mx', 'var_mn', 'angle', 'xE_mean', 'yE_mean', 'varE_mx',
```

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mpixdaq.py

```
class mpixControl:
    """Provides global variables containing Device Information
    as well as Queues, Events and methods to control threads"""
```

```
class miniPIXdaq:
    """Initialize and readout miniPIX device
```

After initialization, the `__call__()` method of this class is executed in an infinite loop, storing data from the device in a ring buffer. The current buffer index is sent to the calling process via a Queue (`dataQ`, an instance of `threading.Queue()`). The loop ends when the flag `*endEvent*` (an instance of `threading.Event()` in class `mpixControl`) is set.

DAQ parameters from class `mpixControl`:

- `ac_count`: number of frames to read successively
- `ac_time`: acquisition (=exposure) time per frame

Queues for communication and synchronization

- `dataQ`: Queue to transfer data
- `cmsQ`: command Queue

Data structure:

- fBuffer: ring buffer with recent frame data

"""

class frameAnalyzer:

"""Analyze frame data and produce a list of cluster objects,

Args: 2d frame data, as obtained from miniPIXdaq.__call__()

Output:

pixel_clusters: a list of tuples of format

(x, y), n_pix, energy, e_mx, (x_mn, y_mn, w, h), (var_mx, var_mn), angle, (xE_m, yE_m), (varE_mx, varE_mn),
with cluster properties

Helper functions to store analysis results are included as static methods

Another static method, cluster_summary() is particularly useful for on-line monitoring of incoming data and provides a summary of the properties of clusters in a pixel frame, returning

- n_clusters: number of multi-pixel clusters
- n_cpixels: number of pixels per cluster
- circularity: circularity per cluster (ranging from 0. for linear, 1. for circular)
- flatness: ratio of maximum variances of pixel and energy distributions in clusters
- cluster_energies: energy per cluster
- single_energies: energies in single pixels

class mpixGraphs:

"""Display of miniPIX frames and histograms for low-rate scenarios, where on-line analysis is possible and animated graphs are meaningful

Animated graph of (overlaid) pixel images, number of clusters per frame and histograms of cluster properties

Args:

- npix: number of pixels per axis (256)
- nover: number of frames to overlay
- unit: unit of energy measurement ("keV" or " μ s ToT")
- circ: circularity of "round" clusters (0. - 1.)
- flat: flatness of energy distribution of pixels in clusters (0. - 1.)
- acq_time: accumulation time per read-out frame
- prescale: prescale factor for frame analysis

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class runDAQ:

"""run miniPIX data acquisition, analysis and real-time graphics and data storage

class to handle:

- command-line arguments
- initialization of miniPIX device of input file

- real-time analysis of data frames
- animated figures to show a live view of incoming data
- event loop controlling data acquisition, data output to file and graphical display

mpixhelpers.py

```
class fileDecoders:
    """Collection of decoders for various input file formats
    supports mPIXdaq .npy and .yaml and Advacam .txt and .clog
    """

class readClusters:
    """Read cluster data written with mPIXdaq, print meta data and statistics

    Default cuts on cluster features are used to classify clusters as  $\alpha$ ,  $\beta$  or  $\gamma$  signatures.

     $\alpha$ : round cluster shape and peaking energy distribution, high ionization per pixel
     $\beta$ : shape is not round, low ionization per pixel, >5 pixels
     $\gamma$ : not ( $\alpha$  or  $\beta$ )

    Methods:

    * __init__(): instantiate class clusterReader and - optionally - set input file name
    * set_cuts(): set cut values
      - small_cut: separate small and large clusters
      - circularity_cut: round topology
      - flatness_cut: flat energy distribution
      - emean_cut: energy loss per pixel (only used if emx not in feature list)
      - emx_cut: maximum pixel energy
      - no_saturation: ignore clusters with saturated pixels
    * parse_args(): read command line arguments if used interactively
    * read_data(): load data in yaml format in pandas data frame
    * set_selection_masks(): define boolean masks to select  $\alpha$ ,  $\beta$  and  $\gamma$ 
    * get_statistics(): count signatures and provide parameters of energy distributions
    * plot(): plot energy distributions of  $\alpha$ ,  $\beta$  and  $\gamma$ 
    * __call__(): execute read_data(), set_selection_masks() and get_statistics()
    """

def pxl2map(pxlist):
    """Create pixel energy map from pixel list

    Args:
    - pxlist: list of pixels [ ..., [px_idx, px_energy], ...]
    - num: int, for numbering figures

    Returns:
    - numpy array
    """

plot_cluster(pxlist, num=0):
    """Plot energy map of pixel cluster

    Args:
    - pxlist: list of pixels [ ..., [px_idx, px_energy], ...]
    - num: int, for numbering figures

    Returns:
    - matplotlib figure
```

```

"""
class shmManager:
    """simple management of shared memory blocks

    class methods:
        - def get_sharedMem(name, size): create or link to shared memory

        do not forget to close() and finally unlink() all requested blocks
        in calling process
    """

```

mplhelpers.py

```

class bhist:
    """one-dimensional histogram for animation, based on bar graph
    supports multiple data classes as stacked histogram

    Args:
        * data: tuple of arrays to be histogrammed
        * bindegges: array of bin edges
        * xlabel: label for x-axis
        * ylabel: label for y axis
        * yscale: "lin" or "log" scale
        * labels: labels for classes
        * colors: colors corresponding to labels
    """

class scatterplot:
    """two-dimensional scatter plot for animation, based on numpy.histogram2d.
    The code supports multiple classes of data and plots a '.' in the corresponding color in every non-
    zero bin of a 2d-histogram

    Args:
        * data: tuple of pairs of coordinates (([x], [y]), ([], []), ...)
        per class to be shown
        * binedges: 2 arrays of bin edges ([bex], [bey])
        * xlabel: label for x-axis
        * ylabel: label for y axis
        * labels: labels for classes
        * colors: colors corresponding to labels
    """

class controlGUI:
    """graphical user interface to control apps via multiprocessing Queue

    Args:
        - cmdQ: a multiprocessing Queue to accept commands
        - appName: name of app to be controlled
        - statQ: mp Queue to show status data
        - confdict: a configuration dictionary for buttons, format {name: [position 0-
5, command]}
    """

```

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```

#!/usr/bin/env python
#
# script run_mPIXdaq.py
# run mpixdaq example with data acquisition, on-line analysis and visualization
# of pixel frames and histogramming

import os
import platform
import sys
import multiprocessing

# on some Linux systems, pypixet requires '.' in LD_LIBRARY_PATH to find C-libraries
# - add current directory to LD-LIBRARY_PATH
# - and restart python script for changes to take effect

modified_path = False
if platform.system() != 'Windows':
    _ldp = os.environ.get("LD_LIBRARY_PATH")
    if _ldp:
        if '://' not in _ldp and _ldp != '.':
            os.environ["LD_LIBRARY_PATH"] = _ldp + ':'
            modified_path = True
    else:
        os.environ['LD_LIBRARY_PATH'] = '.'
        modified_path = True

# restart script in modified environment
if modified_path:
    print(" ! temporarily added '.' to LD_LIBRARY_PATH !")
    try:
        os.execv(sys.argv[0], sys.argv)
    except Exception as e:
        sys.exit('!!! run_mPIXdaq: Failed to Execute under modified environment: ' + str(e))

# get current working directory (before importing minipix libraries)
wd = os.getcwd()
from mpixdaq import mpixdaq # this may change the working directory, depending on system

# finally, start daq in working directory
if __name__ == '__main__': # -----
    if sys.platform.startswith("win"):
        # On Windows calling this function is necessary.
        multiprocessing.freeze_support()

    rD = mpixdaq.runDAQ(wd)
    rD()

```

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```

class materials:
    """Collect properties of target materials and projectiles"""

    dEdx(T, material, projectile):

```

```
"""calculate energy loss in matter ( or "mass stopping power")
for heavy projectiles and electrons
```

Parameters:

```
T: kinetic energy in MeV
Material dictionary
  Z: effective Z
  A: effective A
  I: mean ionization energy (MeV)
z: charge of projectile
m: mass of projectile
```

Returns:

```
mass stopping power in units of MeV cm2/g
(to be multiplied by material density to obtain energy loss per unit length )
```

Formulae:

- Bethe-Bloch relation for heavy projectiles (m_{α} , m_p or m_{μ} $\gg m_e$)
- Notes:
- only collisions considered, no radiation loss (relevant for energy $\gg 1$ MeV)
 - density effect correction (δ) omitted here for simplicity
 - pure Bethe loss decreases below 0.4 MeV and becomes even negative below 0.15 MeV for alpha particles in air; this can be fixed by adding Barkas corrections, which are, however, not implemented here
 - empirical cut-off to avoid negative values of energy loss at small T
- modified Bethe-Bloch equation for electrons, ICRU Report 37 (1984).
- ```
"""
```

```
calc_pixel_energies(E0, px_size=0.0055):
```

```
"""calculate pixel energies for an electron track with energy E0 in silicon
```

```
calc_E_vs_depth(E0, dx, material, projectile):
```

```
"""calculate particle energies after penetrating depth x of material
```

```
plot_dEdx_electron(material, nbins=100, bw=0.05, axis=None):
```

```
"""plot dE/dx * rho for electrons as a function of energy
```

```
plot_beta_pixel_energies(E0=1.0, px_size=0.0055, axis=None):
```

```
"""energy deposits per pixel
```

```
E0: initial energy
```

```
px_size: pixel size in cm
```

```
axis: figure axis; a new one is created if not given
```

```
returns: matplotlib figure
```

```
"""
```

```
plot_dEdx_alpha(material, nbins=200, bw=0.025, axis=None):
```

```
"""dE/dx for alpha particles in material vs. energy
```

```
material: target material
```

```
nbins : number of bins
```

```
bw : bin width
```

```
axis : figure axis; a new one is created if not given
```

```
returns: matplotlib figure
```

```
"""
```

```
plot_alpha_range(material, E0=6.0, dx=0.050, axis=None):
```

```
"""alpha energy after penetration depth in material
```

```
material : target
E0 : initial energy in MeV
dx : step width
axis : figure axis; a new one is created if not given

returns: matplotlib figure
"""
```

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